



Why miss a decibel
of “Tell me more”?

Siemens miniTek™.
One touch.
Endless possibilities.

BestSound™
Technology

The wireless system that unites your audio devices
with your hearing instruments.

Answers for life.

SIEMENS

It all starts with

one touch.

miniTek is a combined streamer and remote control that connects your hearing instruments to multiple audio devices simultaneously.



2 transmitters can be linked to miniTek at once, giving you multiple connections to TVs and stereos.



Bluetooth® mobile phones link automatically to miniTek.



TVs and stereos link to miniTek via the transmitters.



Bluetooth® landline phones are compatible with miniTek. Up to 2 landline and/or mobile phones can work at once.



3rd party Bluetooth® transmitters such as MP3 players and laptops can stream seamlessly to miniTek.



Your **Siemens hearing instruments** communicate wirelessly with miniTek, receiving sound in sync and in stereo.



Integrated Direct Audio Input (DAI) allows use of an **FM transmitter**, and receiver, providing excellent signal-to-noise ratio for learning in a classroom or at home.



Telecoil is built into miniTek, so it can stream from induction loops commonly found in theaters and places of worship – a great advantage if your hearing instruments have no integrated T-coil.

The beauty of one touch.

Both inside and out, miniTek is designed to impress the most discerning tech enthusiast and yet is very simple to operate. With its elegant design and ergonomic buttons, it feels natural from the first touch. And at less than 2 ounces, much lighter than a cell phone, it sits discreetly and comfortably in the palm of your hand.

You've never felt this well-connected.

If you like movies, music, and phone calls, you'll love miniTek™. It unites streamer and remote control functionalities to connect you to all your favorite devices. Enjoy crisp, clear sound from your TV, phone, MP3 player or even computer sent straight to your hearing instruments – completely wirelessly.

Not only does miniTek keep you connected, it can also handle multiple devices at once. Its multipoint system lets you link to 2 phones simultaneously. And it works with 2 transmitters, so you won't have to choose between TVs and your stereo.

Movie night:
82 decibels



Discover a world of smart accessories at the touch of a button.

For over 130 years, our goal has been to ensure that every person enjoys every decibel of life. With solutions that set new standards in the field of hearing instruments. With technologies that ensure efficiency, durability, and respect for the environment. And with accessories that are in tune with your modern lifestyle.

Overview of miniTek & Tek

miniTek

2nd generation Tek remote control and streamer with an ultra-small, high-end design that weighs 2 oz. and is about the size of a matchbox

1 button operation for program selection and volume control

Audio streaming range of 2 feet

Battery life of approximately 5 hours streaming

Offers additional connectivity to 3rd party Bluetooth® wireless technology transmitters, DAI plug for FM receivers, plus Telecoil

Links wearers to up to a total of 7 audio sources, including connection to:

- 2 Bluetooth® enabled devices (mobile and landline phones, laptop)
- 2 Bluetooth® transmitters simultaneously (2 TVs, or your TV and stereo)
- Direct Audio Input (DAI) plug for FM receivers
- T-Coil and a Line-in option (MP3 player)

Tek

Remote control and streamer that provides a clear visual display of battery status, volume, selected programs, and SoundBalance

Large, numbered push buttons for direct selection of programs

Audio streaming range of 3 feet

Battery life of approximately 6 hours streaming

–

Connects wearers to up to a total of 3 devices, including:

- Bluetooth® transmitter (TV)
- Bluetooth® enabled (cell phone)
- Line-in option (MP3 player)



Thai cooking show:
62 decibels

Shawnae Jebbia, Siemens miniTek user and former Miss USA

The brilliance of one touch.

miniTek works with all wireless Siemens hearing instruments to bring you the most natural sound in true stereo. Features like e2e wireless® 2.0 and SoundBrilliance™ ensure clear transmissions for the richest hearing experience.

Remarkable efficiency.

Its rechargeable batteries make miniTek both convenient and environmentally friendly. And compared to similar devices, miniTek offers more functionality without compromising the long battery life of your Siemens hearing instruments – even when you link to 2 devices simultaneously.

Optimal transmission.

To enjoy miniTek's powerful benefits, it should be positioned within 2 feet of your hearing instruments. You can wear it around your neck using a lanyard, or on your clothing with a clip.



Your Audiologist or Hearing Care Provider:

Hearing instruments help many people hear better, but no hearing instrument can solve every hearing problem nor restore normal hearing.

Information in this brochure is subject to change without notice. Copyright © 2011 Siemens Hearing Instruments, Inc. All rights reserved.

The *Bluetooth*® word mark and logos are registered trademarks owned by Bluetooth® SIG, Inc., and any use of such marks by Siemens AG is under license. Other trademarks and trade names are those of their respective owners.

Local Contact Information

Siemens Hearing Instruments, Inc.
P.O. Box 1397
Piscataway, NJ 08855-1397
www.usa.siemens.com/hearing

Global Siemens Headquarters

Siemens AG
Wittelsbacherplatz 2
80333 Muenchen
Germany

Legal Manufacturer

Siemens Audiologische Technik
GmbH
Gebberstrasse 125
91058 Erlangen
Germany

Global Business Unit

Siemens Audiologische Technik
GmbH
Gebberstrasse 125
91058 Erlangen
Germany
Phone: +49 9131 308 0

Global Siemens Healthcare Headquarters

Siemens AG
Healthcare Sector
Henkestrasse 127
91058 Erlangen
Germany
Phone: +49 9131 84 0
www.siemens.com/healthcare

www.usa.siemens.com/minitek